

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen LEPTOSPIRA BIFLEXA

|  |   |
|--|---|
| Best.-Nr. / Order No.                  | 9120  |
| Chargen-Nr. / Lot No.                  | SIH.AU  |
| Verfallsdatum / Expiry Date            | 2027-09   |
| Antigen / Antigen                      | Leptospira biflexa (Stamm Patoc) /<br>Leptospira biflexa /strain Patoc)   |
| Kultursystem / Culture System          | Flüssigkultur / Broth culture   |
| Präparation / Preparation              | mit Ultraschall behandeltes Zellmaterial /<br>sonicated cellular material |
| Inaktivierung / Inactivation           | 0,1 % Thiomersal / 0.1 % thiomersal                                       |
| Konservierung / Preservation           | 0,01 % Thiomersal / 0.01 % thiomersal                                     |
| Lagerung / Storage                     | 2 °C - 8 °C   |
| Rekonstitution / Reconstitution        | 1 ml Aqua dest.   |
| Gebrauchsverdünnung / Working Dilution | 1:8   |



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)