

Analysenzertifikat

Certificate of Analysis

SERION KBR Antigen / CFT Antigen LEPTOSPIRA ICTEROHAEMORRHAGIAE

Best.-Nr. / Order No.	9080
Chargen-Nr. / Lot No.	SGL.AA
Verfallsdatum / Expiry Date	2028-07
Antigen / Antigen	Leptospira interrogans, Serovar icterohaemorrhagiae
Kultursystem / Culture System	Flüssigkultur / Broth culture
Präparation / Preparation	mit Ultraschall behandeltes Zellmaterial / sonicated cellular material
Inaktivierung / Inactivation	0,1 % Thiomersal / 0.1 % thiomersal
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: dialog@virion-serion.de
 Internet: www.virion-serion.de