



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen LISTERIA MONOCYTOGENES

Best.-Nr. / Order No.	1234
Chargen-Nr. / Lot No.	SIF.AP
Verfallsdatum / Expiry Date	2025-09
Antigen / Antigen	Listeria monocytogenes, Serovare I und 4b / Listeria monocytogenes, serovars I and 4b
Kultursystem / Culture System	Plattenkultur / cultivation on agar plates
Präparation / Preparation	beschalltes Ganzzellextrakt / sonicated cell lysate
Inaktivierung / Inactivation	Hitzebehandlung (10 min / 70°C) / heat treatment (10 min / 70°C)
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:30



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 B