

# Analysenzertifikat Certificate of Analysis

## SERION KBR Antigen / CFT Antigen LISTERIA MONOCYTOGENES

Best.-Nr. / Order No.	1234
Chargen-Nr. / Lot No.	SDK.AB
Verfallsdatum / Expiry Date	2029-04
Antigen / Antigen	Listeria monocytogenes, Serovare I und 4b / Listeria monocytogenes, serovars I and 4b
Kultursystem / Culture System	Plattenkultur / cultivation on agar plates
Präparation / Preparation	beschalltes Ganzzellextrakt / sonicated cell lysate
Inaktivierung / Inactivation	Hitzebehandlung (10 min / 70°C) / heat treatment (10 min / 70°C)
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:30



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de