



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Antigen / CFT Antigen MYCOPLASMA PNEUMONIAE

Best.-Nr. / Order No.	1111
Chargen-Nr. / Lot No.	SCF.BH
Verfallsdatum / Expiry Date	2025-03
Antigen / Antigen	Mycoplasma pneumoniae (Stamm FH)/ Mycoplasma pneumoniae (Strain FH)
Kultursystem / Culture System	Flüssigkultur / Broth culture
Präparation / Preparation	Organischer Extrakt aus Mycoplasma pneumoniae, re-extrahiert in eine wässrige Phase / Organic solvent extract of Mycoplasma pneumoniae, re-extracted into an aqueous phase
Inaktivierung / Inactivation	Extraktion der Mycoplasmen mit organischen Lösungsmitteln / Extraction of whole Mycoplasma cells with organic solvents
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:16



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)  
FB-IV-047 B