


Analysenzertifikat Certificate of Analysis



virion\serion

SERION KBR Antigen / CFT Antigen PARAINFLUENZA VIRUS POOL

Best.-Nr. / Order No.	1115
Chargen-Nr. / Lot No.	SDE.CO
Verfallsdatum / Expiry Date	2024-04
Antigen / Antigen	Parainfluenza Type 1 Virus (Stamm/strain Sendai) Parainfluenza Type 2 Virus (Stamm/strain Greer) Parainfluenza Type 3 Virus (Stamm/strain C243)
Kultursystem / Culture System	Parainfluenza Type 1: MDCK-Zellen / MDCK cells Parainfluenza Type 2 und/and 3: Verozellen / Vero cells
Präparation / Preparation	Mischung aus infizierten Zellen und Zellkulturmedium (konzentriert und gereinigt) Mixture of infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Hitzeinaktivierung (56°C) / heat treatment (56°C) Tween/Ether Extraktion / tween/ether extraction
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 A