



# Analysenzertifikat Certificate of Analysis

virion\serion

## SERION KBR Antigen / CFT Antigen PARAINFLUENZA VIRUS POOL

Best.-Nr. / Order No.	1115
Chargen-Nr. / Lot No.	SHF.DG
Verfallsdatum / Expiry Date	2025-08
Antigen / Antigen	Parainfluenza Type 1 Virus (Stamm/strain Sendai) Parainfluenza Type 2 Virus (Stamm/strain Greer) Parainfluenza Type 3 Virus (Stamm/strain C243)
Kultursystem / Culture System	Parainfluenza Type 1: MDCK-Zellen / MDCK cells Parainfluenza Type 2 und/and 3: Verozellen / Vero cells
Präparation / Preparation	Mischung aus infizierten Zellen und Zellkulturmedium (konzentriert und gereinigt) Mixture of infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Hitzeinaktivierung (56°C) / heat treatment (56°C) Tween/Ether Extraktion / tween/ether extraction
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)  
FB-IV-047 B