



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen PARAINFLUENZA VIRUS 3

Best.-Nr. / Order No.	1118
Chargen-Nr. / Lot No.	SKF.BQ
Verfallsdatum / Expiry Date	2025-10
Antigen / Antigen	Parainfluenza Typ 3 Virus /Stamm C243) / Parainfluenza type 3 virus (strain C243)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	Infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) Infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Tween/Ether Behandlung / tween/ether treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 B