



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen PICORNA VIRUS POOL (EV, CAV, CBV)

Best.-Nr. / Order No.	1126
Chargen-Nr. / Lot No.	SGD.BR
Verfallsdatum / Expiry Date	2023-07
Antigen / Antigen	Coxsackie-Viren / Coxsackie viruses A9 , B1, B2, B3, B4, B5, B6 (PB, Conn, Ohio 1, Nancy, JVC, Faulkner, Schmitt) Echo-Viren / Echo viruses 4, 6, 9, 14, 24, 30 (Pesascek, D´Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	LLC-MK ₂ Zellen / cells Vero-Zellen / cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	β-Propiolacton / UV-C
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 A