



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Antigen / CFT Antigen PICORNA VIRUS POOL (EV, CAV, CBV)

Best.-Nr. / Order No.	1126
Chargen-Nr. / Lot No.	SCG.BS
Verfallsdatum / Expiry Date	2026-03
Antigen / Antigen	Coxsackie-Viren / Coxsackie viruses A9, B1, B2, B3, B4, B5, B6 (PB, Conn, Ohio 1, Nancy, JVB, Faulkner, Schmitt) Echo-Viren / Echo viruses 4, 6, 9, 14, 24, 30 (Pesascek, D'Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	LLC-MK <sub>2</sub> Zellen / cells Vero-Zellen / cells
Präparation / Preparation	Infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) Infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	β-Propiolakton / UV-C β-Propiolactone / UV-C
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)  
FB-IV-047 B