



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen ROTAVIRUS

Best.-Nr. / Order No.	1193
Chargen-Nr. / Lot No.	SFC.DQ
Verfallsdatum / Expiry Date	2022-06
Antigen / Antigen	Rotavirus (SA11)
Kultursystem / Culture System	MA-104 Zellen / cells
Präparation / Preparation	Alkalische Extraktion infizierter Zellen (partiell durch Zentrifugation gereinigt) / Alkaline extraction of infected cells (partially purified by centrifugation)
Inaktivierung / Inactivation	0,05 % Triton TX-100 0.05 % Triton TX-100
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de