



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Antigen / CFT Antigen TOXOPLASMA GONDII

Best.-Nr. / Order No.	1331
Chargen-Nr. / Lot No.	SME.AR
Verfallsdatum / Expiry Date	2024-12
Antigen / Antigen	Toxoplasma gondii Stamm / strain RH Toxoplasma gondii Stamm / strain BK
Kultursystem / Culture System	HEL-299
Präparation / Preparation	geerntete und durch Zentrifugation partiell gereinigte Tachyzoiten in Puffer harvested and partially purified tachyzoites in buffer
Inaktivierung / Inactivation	Hitzeinaktivierung / heat inactivation 1 h, 56 °C
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)  
FB-IV-047 B