



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen VARICELLA-ZOSTER VIRUS (VZV)

Best.-Nr. / Order No.	1191
Chargen-Nr. / Lot No.	SFD.BG
Verfallsdatum / Expiry Date	2023-06
Antigen / Antigen	Varizella-Zoster-Virus (Ellen)
Kultursystem / Culture System	HEL299 Zellen / cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) / infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Physikalische Inaktivierung / physical inactivation
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:12

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: dialog@virion-serion.de

Internet: www.virion-serion.de