

Analysenzertifikat

Certificate of Analysis

SERION KBR Antigen / CFT Antigen YERSINIA ENTEROCOLITICA O3

Best.-Nr. / Order No.	1203
Chargen-Nr. / Lot No.	VDI.AB
Verfallsdatum / Expiry Date	2019-10
Antigen / Antigen	Yersinia enterocolitica, serovar O:3
Kultursystem / Culture System	Plattenkultur / cultivation on agar plates
Präparation / Preparation	beschalltes Ganzzellextrakt / sonicated whole cell lysate
Inaktivierung / Inactivation	Hitzebehandlung / Heat treatment (10 min, 70 °C)
Konservierung / Preservation	0,01 % Thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:40



Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: dialog@virion-serion.de
 Internet: www.virion-serion.de