



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Kontrollantigen / CFT Control Antigen CYTOMEGALOVIRUS

Best.-Nr. / Order No.	2130
Chargen-Nr. / Lot No.	SDE.CX
Verfallsdatum / Expiry Date	2024-04
Kultursystem / Culture System	HEL-299
Präparation / Preparation	alkalisch extrahierte Zellen (Kontentriert und gereinigt) cells extracted with alkaline buffer (concentrated and purified)
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: dialog@virion-serion.de

Internet: www.virion-serion.de

serion



multianalyt™