



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Kontrollantigen / CFT Control Antigen COXIELLA BURNETII PHASE 2

Best.-Nr. / Order No.	2123
Chargen-Nr. / Lot No.	SFD.BI
Verfallsdatum / Expiry Date	2023-06
Kultursystem / Culture System	bebrütete Hühnereier incubated hen`s eggs
Präparation / Preparation	Allantoisflüssigkeit des bebrüteten Hühnereis Allantois fluid from incubated hen`s eggs
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10

Hersteller / Manufacturer:

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19

D-97076 Würzburg

Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: dialog@virion-serion.de

Internet: www.virion-serion.de

serion



multianalyt™