



# Analysenzertifikat Certificate of Analysis

virion\serion

## SERION KBR Kontrollantigen / CFT Control Antigen INFLUENZA A/B POOL VIRUS

Best.-Nr. / Order No.	2114
Chargen-Nr. / Lot No.	SFD.BI
Verfallsdatum / Expiry Date	2023-06
Kultursystem / Culture System	bebrütete Hühnereier incubated hen`s eggs
Präparation / Preparation	Allantoisflüssigkeit des bebrüteten Hühnereis  Allantois fluid from incubated hen`s eggs
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)

FB-IV-051 D