

# Institut Virion\Serion GmbH

## Analysenzertifikat / Certificate of Analysis

SERION KBR Kontrollantigen / CFT control antigen  
MASERN VIRUS/MEASLES/ROUGEOLLE  
KONTROLLANTIGEN

Best.-Nr / order no.	2190
Chargen-Nr. / lot no.	SCY.CA
Verfallsdatum / expiry date	2018-03
Zellkulturlinie / Culture system	Verozellen / Vero cells
Material / material	Konzentrierter Zellkulturüberstand und alkalisch extrahierte Zellen concentrated cell culture supernatant and cells extracted with alkaline buffer
Konservierungsmittel / preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / storage	+2°C bis +8°C / +2°C to +8°C
Rekonstitution / reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / working dilution	1:10

16.04.2008



Prüfdatum / Date of Control

Dr. Thomas Schumacher  
Produktspezialist / Product Specialist