



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Kontrollantigen / CFT Control Antigen MUMPS VIRUS/PAROTITIS/OREILLONS

Best.-Nr. / Order No.	2125
Chargen-Nr. / Lot No.	SCC.DL
Verfallsdatum / Expiry Date	2022-03
Kultursystem / Culture System	Verozellen / vero cells
Präparation / Preparation	Konzentrierter Zellkulturüberstand und alkalisch extrahierte Zellen concentrated cell culture supernatant and cells extracted with alkaline buffer
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10

Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)

Internet: [www.virion-serion.de](http://www.virion-serion.de)

serion



multianalyt™