

# Analysenzertifikat

## Certificate of Analysis

**SERION KBR Kontrollantigen / CFT Control Antigen**  
**TOXOPLASMA GONDII**

Best.-Nr. / Order No.	2331
Chargen-Nr. / Lot No.	SMI.DA
Verfallsdatum / Expiry Date	2028-12
Kultursystem / Culture System	HEL-299
Präparation / Preparation	alkalisch extrahierte Zellen (konzentriert und gereinigt)  cells extracted with alkaline buffer (concentrated and purified)
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de