



ut Virion\Serion GmbH · Postfach 92 53 · 97092 Würzburg

Analysenzertifikat Certificate of Analysis

virion\serion
Tel. +49 (0) 931 3045 0
Fax +49 (0) 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

SERION KBR Negativkontrolle / CFT Negative Control PARAINFLUENZA VIRUS POOL

Best.-Nr. / Order No.	4115
Chargen-Nr. / Lot No.	SAB.BG
Verfallsdatum / Expiry Date	2021-01
Material / Material	Fötales Kälberserum / Fetal calf serum
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10
KBR Titer / CFT Titre	< 1:10

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de