

# Analysezertifikat

## Certificate of Analysis

### SERION KBR Negativkontrolle / CFT Negative Control RESP. SYNCYTIAL VIRUS

|  |  |
|--|--|
| Best.-Nr. / Order No.                  | 4124                                   |
| Chargen-Nr. / Lot No.                  | SBI.AD                                 |
| Verfallsdatum / Expiry Date            | 2025-06                                |
| Material / Material                    | Fötales Kälberserum / Fetal calf serum |
| Lagerung / Storage                     | 2 °C - 8 °C                            |
| Rekonstitution / Reconstitution        | 0,1 ml Aqua dest.                      |
| Gebrauchsverdünnung / Working Dilution | 1:10                                   |
| KBR Titer / CFT Titre                  | < 1:10                                 |



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)