

Specifications **SERION easyANALYZE**

To run SERION *easyANALYZE* on your PC, various hardware and software requirements listed below must be met:

- PC with Pentium or AMD processor or higher
- Min. 256 MB main memory
- Min. screen resolution 1024x768
- Serial Port RS-232
- Microsoft® Windows® 7 or 10
- Min. Microsoft® Office 2003, for the use of data export
- Microsoft® .NET Framework 2.0 Redistributable (32 or 64 bit), included in the installation file

For a correct measurement the maintenance and calibration of the photometer according to the manufacturer's instructions is absolutely necessary.

SERION *easyANALYZE* can also be operated with computers with comparable equipment. However, the software is not validated for such PCs.